

Model 8400 Multiplace Double Lock HBOT System



- 84" ID Rectangular Door in Both Locks
- ASME
- Fiber Optics LED Lights
- PVHO-1
- National Board
- NFPA-99- 2015 Compliant
- FDA 510 (k) Cleared
- Medical Air Compressor Package
- Deluge and Hand Line Fire Suppression System
- Seats 9+1 in the Main Lock
- Seats 2+1 in the Entry Lock
- Deluxe Powder Coated Control Console

OPTIONS:

- 6 ATA
- Stretch Version
- ECU (Environmental Control Unit)
- Medical Service Lock
- Patient Handling System

Delivered and Installed in the Continental USA
Lease purchase available



Model 8400 Multiplace Double Lock HBOT System
being loaded for delivery



Tekna, MFG is a world leader in hyperbaric technologies. Purchase of a Tekna HBOT system provides you more than the best hardware in the world; access to exclusive seminars, unrivaled consulting & support, and turnkey solutions to the design and web-presence needs of a provider of hyperbaric services combined make Tekna, MFG the best and only source you need to ensure the success of your venture.