

Model 5100 Multiplace Single Lock HBOT System



- 21st Century Design
- ASME Section VIII
- PVHO-1 (latest addenda)
- NFPA-99-2015
- FDA-510(k)
- 3 ATA Operating Pressure (6 ATA Optional)
- Pressurized with Medical Air
- Deluxe Control Console
- Communications & Entertainment System
- Also available in a Stretch Version at no extra charge

Delivered and Installed in the Continental USA
Lease purchase available



Non-lubricated Medical Air
Package Included



Deluge + Hand-line Fire
Suppression System Included



Tekna MFG is a world leader in hyperbaric technologies. Purchase of a Tekna HBOT system provides you more than the best hardware in the world; access to exclusive seminars, unrivaled consulting & support, and turnkey solutions to the design and web-presence needs of a provider of hyperbaric services combined make Tekna MFG the best and only source you need to ensure the success of your venture.