

Model 7200 DL Double Lock Multiplace HBOT System



- Rectangular door into ML for wheelchair access
- Entry Lock fitted with circular doors
- ASME, PVHO-1, National Board, FDA-510(k) cleared
- Meets or exceeds NFPA-99-2015 requirements
- Delivered in two configurations: 3 ATA and 6 ATA service
- Large acrylic viewports
- Options include a patient handling system, Environmental Control Unit (ECU), and Medical Service Lock.

Delivered and Installed in the Continental USA
Lease purchase available

Note: All TEKNA Multiplace chambers are designed for 6 ATA service but if the buyer only needs 3 ATA service, we do not add the expensive components (gages and valves, etc.) for 6 ATA service. We can upgrade any TEKNA 3 ATA Multiplace system to 6 ATA if required at a later date.



Rectangular Door into ML



Fold-up Seating



Tekna MFG is a world leader in hyperbaric technologies. Purchase of a Tekna HBOT system provides you more than the best hardware in the world; access to exclusive seminars, unrivaled consulting & support, and turnkey solutions to the design and web-presence needs of a provider of hyperbaric services combined make Tekna MFG the best and only source you need to ensure the success of your venture.